



Immunohematology Case Studies

2018 - #8

A Rh(+) patient with Anti-D?

Dr. Fang-Yeh, Chu

*Director, Department of Clinical Pathology & Quality
Management Center, Far Eastern Memorial Hospital*

Email: jacphajacpha@yahoo.com.tw

Clinical History



- Chief Complaint
 - Left thigh pain after accidental fall
- Gender : male
- Age : 9-year-old
- Past History
 - Past medical history : none
 - No significant family history
 - Allergy : no known allergy
 - Transfusion history: none

Clinical History



- Brief History

This 9-year-old boy patient suffered from left thigh pain after accidental fall on 2016-01-04. He was brought to our ER for assessment. The image study revealed left femoral shaft fracture. Orthopedist was consulted and surgical intervention was suggested. Thus, he was admitted to our ward for further evaluation and treatment.

Clinical History



Physical Examination	
General appearance	
Body height :	140cm
Body weight :	35kg
BP :	127 / 59 mmHg
PR :	102 /min
RR :	18 /min
BT :	36.0 °C
Head	
Conjunctiva :	Pink
Sclera :	Anicteric
Pupil :	3.5/3.5
Neck	
Thyroid :	Not enlarged
Lymph node :	No lymphadenopathy

Clinical History



Physical Examination

Thorax

Chest wall :	Symmetric chest expansion
Lung :	Clear breath sound with no rales or wheezing
Breast :	No palpable mass
Heart :	Regular heart beat with no murmur or thrill

Abdomen

Shape :	Floppy
Tenderness :	Nil
Liver :	Impalpable
Spleen :	Impalpable
Scar :	Nil
Palpable mass :	Nil
Kidney :	No knocking pain
Bowel sound :	Normoactive

Clinical History



Physical Examination

Others

Extremities :	Tenderness (+) , deformity and limitation of ROM of left thigh; distal circulation, motor, and sensation: intact
---------------	--

Peripheral pulse

Carotid a.	R : (++) L : (++)
------------	-------------------

Radial a.	R : (++) L : (++)
-----------	-------------------

Dorsalis pedis a.	R : (++) L : (++)
-------------------	-------------------

Femoral a.	R : (++) L : (++)
------------	-------------------

Popliteal a.	R : (++) L : (++)
--------------	-------------------

Posterior tibial a.	R : (++) L : (++)
---------------------	-------------------

Laboratory Data



Test	01-04	01-06	Ref.	Unit
WBC	13.13		3.8~10.4	$10^3/\mu\text{l}$
RBC	4.94		4.1~5.9	$10^6/\mu\text{l}$
HGB	13.9	11.8	13~17	g/dl
PLT	354		140~400	$10^3/\mu\text{l}$
PT	10.0		8.0~12.0	sec
INR	0.98		0.8~1.2	
aPTT	24.7		23.3~39.3	sec

Laboratory Data



Test	01-04	Ref.	Unit
Na	139	135~148	mmol/L
K	3.8	3.5~5.3	mmol/L
Creatinine	0.42	0.7~1.2	mg/dL
AST	26	~37	IU/L
Glucose AC	117	70~100	mg/dL

Clinical Course



- Transfusion Acquisition
 - Preparation for Packed RBC 2 Units at ER on January 4
- Pre-transfusion Test
 - Typing : A+
 - Antibody screening : Positive
 - Autocontrol : 2+
- One unit of Packed RBC was then transfused at OR on January 5 and no transfusion reaction was noted. Post-operative course was uneventful and the patient was discharged

Antibody Identification



Makropanel 16 **LOT** 8000223951 亞東紀念醫院 臨床病理科 血庫組
Makropanel 16-P **LOT** 8000223950 病歷號 _____ 姓名 _____ 血型 **A+** 性別 男 女 輸血史: 無 有(填寫於背面)
Cellbind ID16 **LOT** 8000223959 抗體名稱 _____ P 值 _____ O+ cord cell: _____ Anti-H: _____
Cellbind ID16-P **LOT** n.a. 臨床診斷 _____ DAT: _____ Anti-IgG: _____ Anti-C3d: _____
Column panel 16 **LOT** 8000223955

IVD i 2016-01-21 醫檢師簽章 _____ 操作日期 _____ 主任: _____

Rh-Hr	ID	Rh-Hr													Kell													Duffy		Kidd		Lewis		MNS				Lutheran		Xg		LIAT				
		C	D	E	c	e	C ^m	f	V	K	k	Kp ^a	Kp ^b	Jk ^a	Jk ^b	Js ^a	Js ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Le ^a	Le ^b	P ₁	M	N	S	s	Lu ^a	Lu ^b	Xg ^a	Xg ^b	IS	RT10 ^a	MP	37°C	AHG	CCC								
R ₁ R ₁	1858660	+	+	0	0	+	+	/	/	+	+	0	+	/	+	+	0	+	0	0	0	+	+	+	0	+	+	0	0	+	+	0	+	0	+	0	-	-	2+	-	-	2+				
R ₁ R ₁	177334	+	+	0	0	+	+	/	/	+	+	0	+	/	+	+	0	+	0	0	0	+	+	+	0	+	+	0	0	+	+	0	+	0	+	0	-	-	2+	-	-	2+				
R ₂ R ₂	392234	0	+	+	0	0	0	/	/	0	+	0	+	0	+	0	+	+	+	+	0	0	+	+	+	0	+	+	0	+	+	0	+	+	0	-	-	2+	-	-	2+					
R ₀	301850	0	+	+	+	0	/	/	0	+	0	+	0	+	0	+	+	+	+	0	+	+	+	0	+	+	+	0	+	+	+	+	+	+	+	+	-	-	2+	-	-	2+				
r' ^r	145791	+	0	0	0	+	0	/	/	+	+	0	+	/	+	+	0	+	0	+	0	+	+	+	0	+	+	0	+	+	+	+	+	+	+	+	-	-	-	-	2+					
r' ^r	371033	0	0	+	+	0	0	/	/	0	+	0	+	/	+	+	+	+	0	0	0	+	+	+	0	+	+	0	+	+	+	+	+	+	+	+	-	-	-	-	2+					
rr	300014	0	0	0	+	+	0	/	/	+	0	+	0	+	0	+	+	+	0	0	0	+	+	+	0	+	+	0	+	+	+	+	+	+	+	+	-	-	-	-	2+					
rr	363885	0	0	0	+	+	0	/	/	+	+	0	+	+	+	+	+	0	+	+	+	0	+	+	0	+	+	0	+	+	+	+	+	+	+	+	-	-	-	-	2+					
rr	40200	0	0	0	+	+	0	/	/	+	+	0	+	/	+	+	0	+	+	+	0	+	+	+	0	+	+	0	+	+	+	+	+	+	+	+	-	-	-	-	2+					
rr	302106	0	0	0	+	+	0	/	/	0	+	0	+	0	+	+	0	+	0	0	0	+	+	0	+	+	0	+	+	+	+	+	+	+	+	+	-	-	-	-	2+					
R ₂ R ₁	68923	+	+	+	0	+	0	/	/	0	+	0	+	/	+	+	0	+	0	+	0	+	+	+	0	+	+	0	+	+	+	+	+	+	+	+	-	-	2+	-	-	2+				
R ₂ R ₂	1709905	+	+	+	+	0	0	/	/	0	+	0	+	/	+	0	+	+	+	+	0	+	+	+	0	+	+	+	0	+	+	+	+	+	+	+	-	-	-	-	2+					
rr	1372406	0	0	0	+	+	0	/	/	0	+	+	+	0	+	+	+	0	+	+	+	0	+	+	+	0	+	+	0	+	+	+	+	+	+	-	-	-	-	2+						
r' ^r r	746609	+	0	0	+	+	+	/	/	0	+	0	+	/	+	+	+	+	+	0	+	+	+	0	+	+	+	0	+	+	+	+	+	+	+	+	-	-	-	-	2+					
R ₂ R ₂	15675	0	+	+	+	0	0	/	/	0	+	0	+	0	+	0	+	0	0	0	+	+	+	0	+	+	0	+	+	0	+	+	+	+	+	+	-	-	-	-	2+					
R ₁ R ₁	002410	+	+	0	0	+	0	/	/	+	+	0	+	0	+	+	+	0	+	+	0	0	+	+	+	0	+	+	0	+	+	+	+	+	+	+	-	-	-	-	2+					
Rh-Hr		C	D	E	c	e	C ^m	f	V	K	k	Kp ^a	Kp ^b	Jk ^a	Jk ^b	Js ^a	Js ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Le ^a	Le ^b	P ₁	M	N	S	s	Lu ^a	Lu ^b	Xg ^a	Xg ^b														
抗體篩檢細則組成		Rh-hr													Kell		Duffy		Kidd		Lewis		MNS				P	other		Auto	-	-	2+	-	-	2+										
產品編號: BPLSCE16/P1SCE2		D	C	E	c	e				K	k			Fy ^a	Fy ^b	Jk ^a	Jk ^b	Le ^a	Le ^b	S	s	M	N	Pl	Mi ^a	Dj ^a																				
有效期限: 2016年01月28日		+	0	+	+	0	0	+	0	+	+	0	+	+	+	+	0	+	+	0	+	+	0	+	+	0	0	0	Screen I	-	-	2+	-	-	2+											
		+	+	0	0	+	0	+	0	+	0	+	+	+	0	+	+	0	+	+	0	+	+	0	+	0	0	Screen II	-	-	2+	-	-	2+												
		+	+	0	0	+	0	+	0	+	0	+	+	+	0	+	+	0	+	+	0	+	+	0	+	0	0	Screen III	-	-	2+	-	-	2+												
Selected Cell	Lot No.	PhenoTyping													Result(Titer) ~LIAT~													無法排除之抗體																		
		(A ₂)													(A ₂)													IS	RT10 ^a	MP	37°C	AHG	CCC	1.												
																																		2.												
																																			3.											
																																			4.											
																																			5.											

MP: Manual Polybrene test*
 LIAT: Indirect Antiglobulin Test using Low-ionic-strength saline-albumin enhancement medium

*A rapid manual Polybrene test for detection of red blood cell antibodies have been devised which uses standard laboratory equipment. Red blood cells are incubated with the test sera in a low ionic medium at room temperature for one minute. Polybrene, a quaternary ammonium polymer, is then introduced to cause nonspecific red blood cell aggregation. The test tubes are centrifuged, the cell free supernatant fluid decanted, and the Polybrene effect on the cells is neutralized by adding a dilute sodium citrate-glucose solution.

Interpretations of Identification



- Interpretations
 - Anti-D was identified due to the MP test
 - Negative reactions at the LIAT phase
- Possible situations
 - Anti-D in patient with D variant type?
Patient not transfused
 - Passive anti-D?
No passive anti-D treatment
 - Auto anti-D, possible
 - Auto anti-LW, possible



DTT treatment

History of LW antigen system (ISBT system 016)




- In 1940, Landsteiner and Wiener produced an antibody by injecting rabbits and guinea pigs with RBCs from rhesus monkeys. Thus, they named the antibody “anti-Rhesus”
- A few years later, it was shown that the animal anti-Rhesus were different from the human Rh antibody (i.e., anti-D); thus, the animal anti-Rhesus was renamed anti-LW in honor of Landsteiner and Wiener

Gene, Protein, and Antigens



- The LW glycoprotein, also called *intracellular adhesion molecule 4* (ICAM-4), is a member of the IgSF



Bailly P, Hermand P, Callebaut I, et al. The LW blood group glycoprotein is homologous to intercellular adhesion molecules. *Proc Natl Acad Sci USA* 1994;91:5306–5310.

Niels Lion, David Crettaz, Olivier Rubin, Jean-Daniel Tissot. Stored red blood cells: A changing universe waiting for its map(s). *Journal of Proteomics* 2010;73:374-385.

Gene, Protein, and Antigens




FIGURE 14-1. Model of two proposed membrane complexes containing Band 3 and Rh proteins: 1) containing tetramers of Band 3 and heterotrimers of RhD, RhCE, and RhAG, and linked to the spectrin matrix of the cytoskeleton through Band 3, protein 4.2, and ankyrin; and 2) containing Band 3, RhD, and RhCE, and linked to the spectrin/actin junction through glycophorin C (GPC), p55, and protein 4.1 and through Band 3 and adducin.

Gene, Protein, and Antigens



- The gene (*LW*) encoding the LW glycoprotein consists of three exons distributed over 2.65 kb of DNA on chromosome 19
- There is a phenotypic relationship between LW and the D antigen of the Rh system; D-positive RBCs have stronger expression of LW antigen than D-negative RBCs, and the expression of LW is stronger on cord RBCs than on RBCs from adults

Gene, Protein, and Antigens



- Acquired and often temporary LW-negative phenotypes sometimes occur with production of anti-LW^a or anti-LW^{ab}, a phenomenon that is usually associated with pregnancy or hematologic malignancy
- Such loss of LW antigens is usually associated with the production of LW antibodies
- The strength of LW antigens, decreases from birth until the adult level is reached at about 5 years of age

Gene, Protein, and Antigens



- LW antigens require intramolecular disulfide bonds and the presence of divalent cations, notably Mg^{2+} , for expression
- LW antigens are resistant to treatment of RBCs by ficin, papain, trypsin, α -chymotrypsin (but may be weakened), sialidase, and acid; they are sensitive to treatment of RBCs by pronase and dithiothreitol (DTT)

Antigens of the LW system



Number	Name	Relative frequency	Comments
LW5	LW ^a	High	Antithetical to LW6 (LW ^b); Gln70
LW6	LW ^b	Low	Antithetical to LW5 (LW ^a); Arg70
LW7	LW ^{ab}	High	

Anti-LW Antibodies



- LW antibodies have been reported to be IgM and/or IgG and reactive at room temperature and/or antiglobulin phase
- LW antibodies have been reported to cause mild delayed hemolytic transfusion reaction or mild hemolytic disease of fetus and newborn
- Allo and auto anti LW have been reported

The importance of Anti-LW



- The weak anti-LW may be confused with anti-D auto- or allo-antibodies
- Thus, differentiation between anti-LW and anti-D is important, particularly in women during the fertile period, and in RhD negative-pregnant women, due to the possibility of the need for anti-RhD prophylaxis
- It has been reported that RhD negative blood component was opted for transfusion in the patient with anti-LW

Conclusions



- In this case, the irregular antibody could be possible auto-anti-LW due to non-reactivity with DTT treated RBCS but the antibody was not further characterized and A D- RBC were selected for transfusion

References



- Landsteiner K, Wiener AS. An agglutinable factor in human blood recognized by immune sera for rhesus blood. *Proc Soc Exp Biol Med* 1940;43:223.
- Bailly P, Hermand P, Callebaut I, et al. The LW blood group glycoprotein is homologous to intercellular adhesion molecules. *Proc Natl Acad Sci USA* 1994;91:5306–5310.
- Niels Lion, David Crettaz, Olivier Rubin, Jean-Daniel Tissot. Stored red blood cells: A changing universe waiting for its map(s). *Journal of Proteomics* 2010;73:374-385.
- Mallinson G, Martin PG, Anstee DJ, et al. Identification and partial characterization of the human erythrocyte membrane component(s) that express the antigens of the LW blood-group system. *Biochem J* 1986;234:649–652.
- Bloy C, Hermand P, Blanchard D, et al. Surface orientation and antigen properties of Rh and LW polypeptides of the human erythrocyte membrane. *J Biol Chem* 1990;265:21482–21487
- Reid ME, Oyen R, Marsh WL. Summary of the clinical significance of blood group alloantibodies. *Semin Hematol* 2000;37:197–216.
- Daniels G. *Human Blood Groups*. 2nd ed. Oxford, England: Blackwell Science; 2002.
- Miola MP, Cervo SV, Fachini RM, Ricci Júnior O. Do not confuse anti-LW autoantibodies with anti-D. *Rev Bras Hematol Hemoter*. 2013;35:198.
- Roback J, Combs MR, Grossman B, et al. *Technical Manual*. 16th ed. Washington, DC: *American Association of Blood Banks*; 2008.